

# RESTART

ULTRASOUND TREATMENT

## Sports Device

**RESTART** device for ultrasound treatment is ideal for all muscle, tendon and ligament treatments.

It is also suggested using the device for the following aesthetic treatments: cellulite, draining massage and microlifting.

The device can be used by the user himself or by a professional, for example an athletic trainer.



**DIMENSIONS**  
204 x 63 x 58 mm

POWERED BY  
**I-TECH**  
MEDICAL DIVISION

## EQUIPMENT

- Portable device
- Power supply
- Battery Ni-Mh AAA850mAh 4.8V
- Ultrasound gel
- User manual
- Carriage and storage bag

## TECHNICAL FEATURES

- Ultrasound treatment device
- 1MHz ultrasound, maximum output power 6.4W (with power supply) and 3.2W (with battery)
- Maximum intensity 1.6W/cm<sup>2</sup> (with power supply) and 0.8W/cm<sup>2</sup> (with battery)
- Collimated ultrasound beam for an optimal penetration through tissues
- 3 ultrasound beam intensity levels: L (low), M (medium) e H (high)
- Medical power supply and battery mode
- Direct contact device (not suitable for water immersion)

SK459-02

## TREATMENTS

Muscle

Tendon/Ligament

Aesthetic

**WEPERE**

IACER srl - via Enzo Ferrari 2 - 30037 Scorzè (VE) - Italia



[www.wepere.com](http://www.wepere.com)