

RESTART

ULTRASOUND TREATMENT

Sports Device

RESTART device for ultrasound treatment is ideal for all muscle, tendon and ligament treatments.

It is also suggested using the device for the following aesthetic treatments: cellulite, draining massage and microlifting.

The device can be used by the user himself or by a professional, for example an athletic trainer.



CND Z120601
REP 5722/R
MEDICAL DEVICE CLASS IIa
CE0068-certified medical device

DIMENSIONS
204 x 63 x 58 mm

POWERED BY
I-TECH
MEDICAL DIVISION

EQUIPMENT

- Portable device
- Power supply
- Battery Ni-Mh AAA850mAh 4.8V
- Ultrasound gel
- User manual
- Carriage and storage bag

TECHNICAL FEATURES

- Ultrasound treatment device
- 13 programs
- 1MHz ultrasound, maximum output power 6.4W (with power supply) and 3.2W (with battery)
- Maximum intensity 1.6W/cm² (with power supply) and 0.8W/cm² (with battery)
- Collimated ultrasound beam for an optimal penetration through tissues
- 3 ultrasound beam intensity levels: L (low), M (medium) e H (high)
- Medical power supply and battery mode
- Direct contact device (not suitable for water immersion)

SK459-01

TREATMENTS

Muscle

Tendon/Ligament

Aesthetic

WEPERE

IACER srl - via Enzo Ferrari 2 - 30037 Scorzè (VE) - Italia



www.wepere.com