



Pressotherapy acts as a lymphatic drainage massage that promotes proper blood and lymphatic circulation and the drainage of liquids and toxins.

Ideal for reducing recovery time after physical activity and improving performance.

ARYA Limited Edition Edition is dedicated to the most ambitious sportspeople who want to stand out.

A hand-crafted and completely **Made in Italy**.



INDICATED AND SUGGESTED TO

- Muscle recovery
- Drainage of liquids
- Improvement of the blood and lymphatic circulation

TECHNICAL FEATURES

Power supply	230AC, 50 Hz
Total consumption	100mA MAX
Therapy time.....	0 ÷ 30 minutes
Pressure	200 mmHg (±20%)
Weight.....	2 kg

WHICH ARE THE INTERMITTENT PNEUMATIC COMPRESSION EFFECTS?

Intermittent pneumatic compression is a useful method for improving a correct venous circulation, reducing muscle tension caused by stress or chronic and acute pain.

Thanks to its pump action, it promotes venous return, increasing the blood supply of the tissues and their consequent correct physiological renewal.

External compression also allows excess interstitial fluids to return more easily to the circulatory system, so that they can be properly and quickly eliminated.

EQUIPMENT

2 boots, 2 connection tubes with 4 terminals, 1 connector, 1 splitter.

Boots also available in size M or L (ARYA Limited Edition LEG2 M or ARYA Limited Edition LEG2 L).

Also available as optional item is the KIT BEAUTY (abdominal band, connection tube with 4 terminals, splitter) for simultaneous treatment of legs and abdomen.

SK620-00